



# DNAJB1 Antibody Pair

<b>Product Code</b>	CSB-EAP03934
<b>Species Reactivity</b>	Human,Rat,Hamster,Mouse,Guinea pig,Bovine,Horse,Pig,Dog,Rabbit,Goat,Sheep
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	DnaJ protein homolog 1,Heat shock 40 kDa protein 1,HSP40,Heat shock protein 40,Human DnaJ protein 1,hDj-1,DNAJB1,DNAJ1, HDJ1, HSPF1
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	DNAJB1
<b>Protein Names</b>	DnaJ homolog subfamily B member 1
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP03934C Detection: CSB-EAP03934D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.