



RETN Antibody Pair

Product Code	CSB-EAP07079
Species Reactivity	Human,Mouse,Guinea pig,Horse,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	Adipose tissue-specific secretory factor,ADSF,C/EBP-epsilon-regulated myeloid-specific secreted cysteine-rich protein,Cysteine-rich secreted protein A12-alpha-like 2,Cysteine-rich secreted protein FIZZ3,UNQ407/PRO1199,FIZZ3, HXCP1, RSTN,RETN
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RETN
Protein Names	Resistin
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07079C Detection: CSB-EAP07079D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.