



PSME2 Antibody Pair

Product Code	CSB-EAP01134
Species Reactivity	Human,Mouse,Hamster,Bovine,Pig,Dog,Chicken,Rabbit,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	PSME2,11S regulator complex subunit beta,REG-beta Activator of multicatalytic protease subunit 2, Proteasome activator 28 subunit beta,PA28b Short name=PA28beta
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	PSME2
Protein Names	Proteasome activator complex subunit 2
Notes	We recommend using the capture antibody at a concentration of 1.6ug/ml and the detection antibody at a concentration of 0.34ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01134C Detection: CSB-EAP01134D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.