



PSMB2 Antibody Pair

Product Code	CSB-EAP01385
Species Reactivity	Human,Mouse,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	Macropain subunit C7-I,Multicatalytic endopeptidase complex subunit C7-I,Proteasome component C7-I
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	PSMB2
Protein Names	Proteasome subunit beta type-2
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.125ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01385C Detection: CSB-EAP01385D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4