



BCAP31 Antibody Pair

Product Code	CSB-EAP03425
Species Reactivity	Human, Mouse, Guinea pig, Bovine, Dog, Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	BCR-associated protein 31, Bap31, 6C6-AG tumor-associated antigen, Protein CDM, p28, BCAP31, BAP31, DXS1357E
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	BCAP31
Protein Names	B-cell receptor-associated protein 31
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.125ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP03425C Detection: CSB-EAP03425D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.