



RPS16 Antibody Pair

Product Code	CSB-EAP02465
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Dog,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	RPS16
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RPS16
Protein Names	40S ribosomal protein S16
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP02465C Detection: CSB-EAP02465D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.