



MRPL9 Antibody Pair

Product Code	CSB-EAP04864
Species Reactivity	Human, Mouse, Guinea pig, Bovine, Pig, Dog, Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	L9mt, MRP-L9, MRPL9
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	MRPL9
Protein Names	39S ribosomal protein L9, mitochondrial
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP04864C Detection: CSB-EAP04864D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.