



NDUFV2 Antibody Pair

Product Code	CSB-EAP05385
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Horse,Pig,Dog,Rabbit,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	NADH-ubiquinone oxidoreductase 24 kDa subunit,NDUFV2
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	NDUFV2
Protein Names	NADH dehydrogenase [ubiquinone] flavoprotein 2, mitochondrial
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP05385C Detection: CSB-EAP05385D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4