



RPS14 Antibody Pair

Product Code	CSB-EAP04255
Species Reactivity	Human, Mouse, Guinea pig, Bovine, Pig, Dog, Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	RPS14 PRO2640
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RPS14
Protein Names	40S ribosomal protein S14
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP04255C Detection: CSB-EAP04255D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.