



HNRPH1 Antibody Pair

Product Code	CSB-EAP04264
Species Reactivity	Human,Mouse,Guinea pig,Bovine,Horse,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	hnRNP H,HNRNPH1,HNRPH, HNRPH1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	HNRNPH1
Protein Names	Heterogeneous nuclear ribonucleoprotein H
Notes	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP04264C Detection: CSB-EAP04264D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4