



# MFGE8 Antibody Pair

<b>Product Code</b>	CSB-EAP07974
<b>Species Reactivity</b>	Human, Guinea pig, Bovine, Horse, Pig, Dog, Goat
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	Breast epithelial antigen BA46 HMFG MFGM Milk fat globule-EGF factor 8 Short name=MFG-E8 SED1 Cleaved into the following 2 chains: Lactadherin short form Medin
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	MFGE8
<b>Protein Names</b>	Lactadherin
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP07974C Detection: CSB-EAP07974D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.