



MFGE8 Antibody Pair

Product Code	CSB-EAP07974
Species Reactivity	Human, Guinea pig, Bovine, Horse, Pig, Dog, Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Breast epithelial antigen BA46 HMFG MFGM Milk fat globule-EGF factor 8 Short name=MFG-E8 SED1 Cleaved into the following 2 chains: Lactadherin short form Medin
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	MFGE8
Protein Names	Lactadherin
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07974C Detection: CSB-EAP07974D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4