



BTF3 Antibody Pair

Product Code	CSB-EAP01714
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Horse,Pig,Dog,Chicken,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	Nascent polypeptide-associated complex subunit beta,NAC-beta,RNA polymerase B transcription factor 3,BTF3,NACB,OK/SW-cl.8
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	BTF3
Protein Names	Transcription factor BTF3
Notes	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01714C Detection: CSB-EAP01714D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.