



# VEGFC Antibody Pair

<b>Product Code</b>	CSB-EAP07545
<b>Species Reactivity</b>	Human,Rat,Guinea pig,Bovine,Horse,Pig,Dog,Goat,Sheep,Duck
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	VEGF-C,Flt4 ligand,Flt4-L,Vascular endothelial growth factor-related protein,VRP,VEGFC
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Gene Names</b>	VEGFC
<b>Protein Names</b>	Vascular endothelial growth factor C
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP07545C Detection: CSB-EAP07545D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4