



CSNK1E Antibody Pair

Product Code	CSB-EAP01315
Species Reactivity	Human,Chicken
Tested Applications	S-ELISA
Form	Liquid
Alias	CSNK1E,CKI-epsilon,CKIe
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	CSNK1E
Protein Names	Casein kinase I isoform epsilon
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01315C Detection: CSB-EAP01315D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.