



PSMB4 Antibody Pair

Product Code	CSB-EAP01374
Species Reactivity	Human, Guinea pig, Pig, Chicken, Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	26 kDa prosomal protein HsBPROS26 PROS-26 Macropain beta chain Multicatalytic endopeptidase complex beta chain Proteasome chain 3 HsN3
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	PSMB4
Protein Names	Proteasome subunit beta type-4
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01374C Detection: CSB-EAP01374D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.