



SNRPG Antibody Pair

Product Code	CSB-EAP01535
Species Reactivity	Human,Pig,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	snRNP-G?Sm protein G?Sm-G?SmG
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	SNRPG
Protein Names	Small nuclear ribonucleoprotein G
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP01535C Detection: CSB-EAP01535D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.