



SAT1 Antibody Pair

Product Code	CSB-EAP02754
Species Reactivity	Human, Guinea pig, Bovine, Pig, Dog, Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	Polyamine N-acetyltransferase 1, Putrescine acetyltransferase, Spermidine/spermine N(1)-acetyltransferase 1, SSAT, SSAT-1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	SAT1
Protein Names	Diamine acetyltransferase 1
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.1ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP02754C Detection: CSB-EAP02754D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.