



NME1 Antibody Pair

Product Code	CSB-EAP02775
Species Reactivity	Human,Guinea pig,Bovine,Pig,Dog,Chicken
Tested Applications	S-ELISA
Form	Liquid
Alias	NDK A,NDP kinase A,A-activated DNase,GAAD,Metastasis inhibition factor nm23,Tumor metastatic process-associated protein,nm23-H1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	NME1
Protein Names	Nucleoside diphosphate kinase A
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP02775C Detection: CSB-EAP02775D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.