



RPS6KB2 Antibody Pair

Product Code	CSB-EAP02905
Species Reactivity	Human,Chicken
Tested Applications	S-ELISA
Form	Liquid
Alias	S6K-beta-2,S6K2,70 kDa ribosomal protein S6 kinase 2,S6 kinase-related kinase,Serine/threonine-protein kinase 14B,p70 ribosomal S6 kinase beta,S6K-beta,p70 S6 kinase beta,p70 S6K-beta,RPS6KB2,STK14B
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RPS6KB2
Protein Names	Ribosomal protein S6 kinase beta-2
Notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP02905C Detection: CSB-EAP02905D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.