



PRMT1 Antibody Pair

Product Code	CSB-EAP03025
Species Reactivity	Human
Tested Applications	S-ELISA
Form	Liquid
Alias	Histone-arginine N-methyltransferase PRMT1 Interferon receptor 1-bound protein 4
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	PRMT1
Protein Names	Protein arginine N-methyltransferase 1
Notes	We recommend using the capture antibody at a concentration of 2ug/ml and the detection antibody at a concentration of 0.2ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP03025C Detection: CSB-EAP03025D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.