



EXTL3 Antibody Pair

Product Code	CSB-EAP03075
Species Reactivity	Human, Mouse, Guinea pig, Bovine, Pig, Dog, Chicken, Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	EXT-related protein 1 Glucuronyl-galactosyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase Hereditary multiple exostoses gene isolog Multiple exostosis-like protein 3 Putative tumor suppressor protein EXTL3
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	EXTL3
Protein Names	Exostosin-like 3
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP03075C Detection: CSB-EAP03075D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.