



RPS10 Antibody Pair

Product Code	CSB-EAP02565
Species Reactivity	Human,Guinea pig,Pig,Dog,Chicken
Tested Applications	S-ELISA
Form	Liquid
Alias	RPS10
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	RPS10
Protein Names	40S ribosomal protein S10
Notes	We recommend using the capture antibody at a concentration of 1ug/ml and the detection antibody at a concentration of 0.125ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP02565C Detection: CSB-EAP02565D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.