



B2M Antibody Pair

Product Code	CSB-EAP002486
Species Reactivity	Human,Rat,Mouse,Guinea pig,Bovine,Pig,Dog,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	B2M CDABP0092, HDCMA22P
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	B2M
Protein Names	Beta-2-microglobulin
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP002486C Detection: CSB-EAP002486D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4