



ARHGEF18 Antibody Pair

Product Code	CSB-EAP05275
Species Reactivity	Human,Mouse,Pig
Tested Applications	S-ELISA
Form	Liquid
Alias	114 kDa Rho-specific guanine nucleotide exchange factor,p114-Rho-GEF,p114RhoGEF,Septin-associated RhoGEF,SA-RhoGEF,ARHGEF18,KIAA0521
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	ARHGEF18
Protein Names	Rho guanine nucleotide exchange factor 18
Notes	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP05275C Detection: CSB-EAP05275D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.