



CSRP2 Antibody Pair

Product Code	CSB-EAP06937
Species Reactivity	Human,Mouse,Guinea pig,Horse,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	Cysteine-rich protein 1,CRP,CRP1,Epididymis luminal protein 141,HEL-141,CSRP1,CSRP, CYRP
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	CSRP1
Protein Names	Cysteine and glycine-rich protein 1
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP06937C Detection: CSB-EAP06937D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.