



FOXP3 Antibody Pair

Product Code	CSB-EAP15229
Species Reactivity	Human,Rat,Horse,Pig,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	JM2,IPEX,FOXP3,Scurfin
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	FOXP3
Protein Names	Forkhead box protein P3
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP15229C Detection: CSB-EAP15229D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4