



# APP Antibody Pair

Product Code	CSB-EAP07174
Species Reactivity	Human,Rat,Mouse,Bovine,Pig,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	ABPP APPI APP Alzheimer disease amyloid protein Cerebral vascular amyloid peptide CVAP PreA4 Protease nexin-II PN-II
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	APP
Protein Names	Amyloid beta A4
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07174C Detection: CSB-EAP07174D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4