



TTR Antibody Pair

Product Code	CSB-EAP12077
Species Reactivity	Human,Rat,Hamster,Mouse,Guinea pig,Bovine,Horse,Pig,Dog,Chicken,Rabbit,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	ATTR,Prealbumin,TBPA,TTR,PALB
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	TTR
Protein Names	Transthyretin
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP12077C Detection: CSB-EAP12077D(HRP) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: Preservative: 0.03% Proclin 300Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.