



MMP1 Antibody Pair

Product Code	CSB-EAP07009
Species Reactivity	Human,Rat,Mouse,Guinea pig,Horse,Pig,Dog,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Fibroblast collagenase,Matrix metalloproteinase-1,MMP-1
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	MMP1
Protein Names	Interstitial collagenase
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP07009C Detection: CSB-EAP07009D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.