



IL18 Antibody Pair

Product Code	CSB-EAP08377
Species Reactivity	Human,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	IGIF, IL1F4,IL18,Iboctadekin,Interferon gamma-inducing factor,IFN-gamma-inducing factor,Interleukin-1 gamma,IL-1 gamma
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	IL18
Protein Names	Interleukin-18
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP08377C Detection: CSB-EAP08377D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.