



PPBP Antibody Pair

Product Code	CSB-EAP09014
Species Reactivity	Human, Mouse, Guinea pig, Bovine, Pig
Tested Applications	S-ELISA
Form	Liquid
Alias	C-X-C motif chemokine 7, LDGF, MDGF, Small-inducible cytokine B7, CTAP3, CXCL7, SCYB7, TGB1, THBGB1, PPBP
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	PPBP
Protein Names	Platelet basic =
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.34ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP09014C Detection: CSB-EAP09014D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.