



Cxcl1 Antibody Pair

Product Code	CSB-EAP09217
Species Reactivity	Human,Rat,Mouse,Guinea pig,Pig,Duck
Tested Applications	S-ELISA
Form	Liquid
Alias	C-X-C motif chemokine 1,Platelet-derived growth factor-inducible protein KC,Secretory protein N51,Hematopoietic synergistic factor,HSF,KC-T,Cxcl1,Gro, Gro1, Mgsa, Scyb1
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Gene Names	Cxcl1
Protein Names	Growth-regulated alpha
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP09217C Detection: CSB-EAP09217D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4