



Cxcl5 Antibody Pair

Product Code	CSB-EAP09257
Species Reactivity	Rat,Mouse,Pig
Tested Applications	S-ELISA
Form	Liquid
Alias	Cytokine LIX Small-inducible cytokine B5
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Target Names	Cxcl5
Protein Names	C-X-C motif chemokine 5
Notes	We recommend using the capture antibody at a concentration of 0.2ug/ml and the detection antibody at a concentration of 0.2ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP09257C Detection: CSB-EAP09257D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.