



BDNF Antibody Pair

Product Code	CSB-EAP05775
Species Reactivity	Human, Guinea pig, Bovine, Pig, Dog, Rabbit, Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Abrineurin, BDNF
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	BDNF
Protein Names	Brain-derived neurotrophic factor
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP05775C Detection: CSB-EAP05775D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4