



Rb1 Antibody Pair

Product Code	CSB-EAP15529
Species Reactivity	Rat,Pig,Dog,Rabbit,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	pRb,Rb,pp105,Rb1,Rb-1
Product Type	Antibody Pairs
Immunogen Species	Rattus norvegicus (Rat)
Gene Names	Rb1
Protein Names	Retinoblastoma-associated protein
Notes	We recommend using the capture antibody at a concentration of 0.3ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP15529C Detection: CSB-EAP15529D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4