



Egln3 Antibody Pair

Product Code	CSB-EAP14797
Species Reactivity	Rat,Mouse,Horse,Pig,Chicken,Rabbit,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Hypoxia-inducible factor prolyl hydroxylase 3,Prolyl hydroxylase domain-containing protein 3,SM-20,Egln3
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Target Names	Egln3
Protein Names	Egl nine homolog 3
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP14797C Detection: CSB-EAP14797D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.