



# Hp Antibody Pair

<b>Product Code</b>	CSB-EAP17069
<b>Species Reactivity</b>	Rat,Guinea pig,Bovine,Horse,Pig,Dog,Chicken,Rabbit,Goat,Sheep
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	Liver regeneration-related protein LRRG173, Zonulin,Hp,Ba1-647
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Target Names</b>	Hp
<b>Protein Names</b>	Haptoglobin
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 1ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP17069C Detection: CSB-EAP17069D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.