



MT1E Antibody Pair

| | |
|----------------------------|---|
| Product Code | CSB-EAP17364 |
| Species Reactivity | Human,Bovine,Horse,Pig,Dog,Rabbit,Goat |
| Tested Applications | S-ELISA |
| Form | Liquid |
| Alias | Metallothionein-IE,MT1E |
| Product Type | Antibody Pairs |
| Immunogen Species | Homo sapiens (Human) |
| Gene Names | MT1E |
| Protein Names | Metallothionein-1E |
| Notes | We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher. |
| Host | Capture: Rabbit Detection: Rabbit |
| Components | Capture: CSB-EAP17364C Detection: CSB-EAP17364D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol. |
| Storage-Buffer | Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4 |