



Muc5ac Antibody Pair

Product Code	CSB-EAP17229
Species Reactivity	Mouse?Guinea pig?Pig?Dog?Rabbit?Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	Muc5ac
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Gene Names	Muc5ac
Protein Names	Protein Muc5ac
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP17229C Detection: CSB-EAP17229D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4