



ANG1 Antibody Pair

Product Code	CSB-EAP14679
Species Reactivity	Bovine,Horse,Pig,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	ANG1,ANG
Product Type	Antibody Pairs
Immunogen Species	Bos taurus (Bovine)
Target Names	ANG1
Protein Names	Angiogenin-1
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP14679C Detection: CSB-EAP14679D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.