



IGF1R Antibody Pair

Product Code	CSB-EAP08399
Species Reactivity	Human,Bovine,Horse,Pig,Dog,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Insulin-like growth factor I receptor,CD_antigen=CD221,IGF1R
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	IGF1
Protein Names	Insulin-like growth factor 1 receptor
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP08399C Detection: CSB-EAP08399D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4