



IL1b Antibody Pair

Product Code	CSB-EAP17347
Immunogen	Rabbit
Species Reactivity	Human?Mouse?Bovine?Pig?Goat
Tested Applications	S-ELISA
Form	Liquid
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Protein Names	Interleukin-1 beta
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml. Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP17347C Detection: CSB-EAP17347D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4