



NEU3 Antibody Pair

Product Code	CSB-EAP17267
Species Reactivity	Human,Pig
Tested Applications	S-ELISA
Form	Liquid
Alias	Ganglioside sialidase,Membrane sialidase,N-acetyl-alpha-neuraminidase 3,NEU3
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	NEU3
Protein Names	Sialidase-3
Notes	We recommend using the capture antibody at a concentration of 0.1ug/ml and the detection antibody at a concentration of 0.125ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP17267C Detection: CSB-EAP17267D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4