



# LAMB2 Antibody Pair

Product Code	CSB-EAP14279
Species Reactivity	Human,Bovine,Horse,Pig,Dog,Chicken,Goat
Tested Applications	S-ELISA
Form	Liquid
Alias	Laminin B1s chain,Laminin-11 subunit beta,Laminin-14 subunit beta,Laminin-15 subunit beta.LAMB2 ,LAMS
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Gene Names	LAMB2
Protein Names	Laminin subunit beta-2
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.5ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP14279C Detection: CSB-EAP14279D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4