



Alpl Antibody Pair

Product Code	CSB-EAP17039
Species Reactivity	Mouse,Bovine,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	Alkaline phosphatase 2,Alpl,Akp-2, Akp2
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Gene Names	Alpl
Protein Names	Alkaline phosphatase, tissue-nonspecific isozyme
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP17039C Detection: CSB-EAP17039D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4