



# PTH Antibody Pair

<b>Product Code</b>	CSB-EAP09777
<b>Species Reactivity</b>	Human,Mouse,Guinea pig,Bovine,Pig,Dog,Goat
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	Parathormone,Parathyrin,PTH
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Protein Names</b>	Parathyroid hormone
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.2ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP09777C Detection: CSB-EAP09777D
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4