



# ASGR1 Antibody Pair

<b>Product Code</b>	CSB-EAP17469
<b>Species Reactivity</b>	Human,Bovine,Pig,Dog
<b>Tested Applications</b>	S-ELISA
<b>Form</b>	Liquid
<b>Alias</b>	C-type lectin domain family 4 member H1,Hepatic lectin H1, ASGR1,CLEC4H1
<b>Product Type</b>	Antibody Pairs
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Target Names</b>	ASGR1
<b>Protein Names</b>	Asialoglycoprotein receptor 1
<b>Notes</b>	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
<b>Host</b>	Capture: Rabbit Detection: Rabbit
<b>Components</b>	Capture: CSB-EAP17469C Detection: CSB-EAP17469D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
<b>Storage-Buffer</b>	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
<b>Usage</b>	For Research Use Only. Not for use in diagnostic or therapeutic procedures.