



HDAC7 Antibody Pair

Product Code	CSB-EAP17899
Species Reactivity	Human,Bovine,Pig,Dog,Goat,Sheep
Tested Applications	S-ELISA
Form	Liquid
Alias	Histone deacetylase 7A,HDAC7,HDAC7A
Product Type	Antibody Pairs
Immunogen Species	Homo sapiens (Human)
Target Names	HDAC7
Protein Names	Histone deacetylase 7
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.34ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP17899C Detection: CSB-EAP17899D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.