



IL18bp Antibody Pair

Product Code	CSB-EAP16769
Species Reactivity	Mouse,Pig,Dog
Tested Applications	S-ELISA
Form	Liquid
Alias	Tadekinig-alfa,IL18BP
Product Type	Antibody Pairs
Immunogen Species	Mus musculus (Mouse)
Target Names	IL18bp
Protein Names	Interleukin-18-binding
Notes	We recommend using the capture antibody at a concentration of 0.5ug/ml and the detection antibody at a concentration of 0.25ug/ml.Optimal dilutions should be determined experimentally by the researcher.
Host	Capture: Rabbit Detection: Rabbit
Components	Capture: CSB-EAP16769C Detection: CSB-EAP16769D(Biotin) Reagents are sufficient for at least 5 x 96 well plates using recommended protocol.
Storage-Buffer	Capture: 50% Glycerol, 0.01M PBS, PH 7.4 Detection: 50% Glycerol, 0.01M PBS, PH 7.4
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.